MMR

HIGH PRESSURE REGULATORS

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
KÄYTTÖOHJE

1. FOREWORD

GCE Medical Regulators are medical devices classified as class IIb according to the Medical Device Directive 93/42/EEC. Their Compliance with essential requirements of 93/42/EEC Medical Device Directive is based upon ISO 10524-1 standard.

2. INTENDED USE

GCE Medical Regulators are designed for use with high-pressure medical gas cylinders equipped with a medical cylinder shut-off valve and in pipeline systems. They regulate pressure and flow of medical gases to the patient. They are intended for the administration of the following medical gases in the treatment, management, diagnostic evaluation and care of the patient.

Basic variants (divided by gas type):

- oxygen;
- xenon;
- nitrous oxide;
- specified mixtures of the gases listed;
- air for breathing;
- air to power surgical tools;
- helium;
- nitrogen to power surgical tools;
- carbon dioxide

3. OPERATIONAL, TRANSPORT AND STORAGE SAFETY REQUIREMENTS

Keep the product and its associated equipment from:

- heat sources (fire, cigarettes,...),
- flammable materials,
- oil or grease (especially be careful in use of hand cream),
- water,
- dust.

The product and its associated equipments must be prevented from falling over.

Always maintain oxygen cleanliness standards.

Use the product and its associated equipment in well ventilated areas only.

Before initial use, the product must be kept in its original packaging. GCE recommends use of the original packaging (including internal sealing bag and caps) if the product is withdrawn from operation (for transport, storage).

Statutory laws, rules, and regulations for medical gases, occupational safety, and environmental protection must be observed.


OPERATING CONDITIONS		STORAGE AND TRANSPORT CONDITIONS	
	-20/+60 °C		-30/+60 °C
	10/100 %		10/100 %
	600/1200 mbar		600/1200 mbar

4. PERSONNEL INSTRUCTIONS

The Medical Devices Directive 93/42/EEC states that product provider must ensure that all personnel handling the product are provided with the operating instructions & performance data.

 **Do not use the product without properly familiarization of the product and its safe operation as defined in this Instruction for use. Ensure user is aware of particular information and knowledge required for the gas in use.**

5. PRODUCT DESCRIPTION

A - Inlet stem

The regulator is fitted to the medical cylinder valve and to the pipeline system by means of an inlet stem. The stem's nut can be female or male threaded. The inlet stem contains a filter.

B - Inlet pressure gauge

The pressure gauge is designed to indicate cylinder contents (the medical cylinder shut-off valve must be set to the ON position to allow pressure measuring).

C - Outlet pressure gauge

The pressure gauge is designed to indicate outlet pressure in the valve and in the pipeline system.

D - Pressure outlet

The regulator may be fitted with a pressure outlet. The pressure outlet is a direct outlet from the low-pressure chamber. One type of pressure outlet is used: The pressure outlet is fitted with a threaded connector. Regulators with this type of pressure outlet can only be used as an integral part of a unit of medical equipment (e.g. emergency ventilator, anaesthesia device, etc.) or in a pipeline system.

E - Pressure relief valve

Function of the pressure relief valve is to relieve excessive overpressure in the valve.

6. OPERATIONS

6.1. BEFORE USE

Visual Inspection before use

- Check whether there is visible external damage to the regulator or to the gas cylinder (including product labels and marking). If they show signs of external damages, remove them from service and identify their status.
- Check visually whether the regulator or the medical gas cylinder is contaminated; if needed, clean the regulator in accordance with the cleaning procedure detailed hereinafter, (in case of cylinder contamination, refer to the gas cylinder manufacturer cleaning recommendations).
- Check whether the relevant product service term or the total life time of the GCE product and of the gas cylinder have not been exceeded (refer to date coding system of GCE or of the owner). If the service term or the life time has been exceeded, remove the regulator (or the gas cylinder) from service & identify its status appropriately.
- Ensure that the product inlet stem is compatible with the medical cylinder valve (gas/ thread type).
- Check the presence & integrity of the inlet stem seals / correct size of seal.

-  **Remove the protection cap from the inlet and/or flow outlet. Keep the caps in a safe place for possible reuse during transport or storage. The product is intended for use only with the gas specified on its label. Never try to use it for another gas!**

Fitting to the medical cylinder valve

- Secure the gas cylinder's position. signs of external damages, remove them from service and identify their status.

Screw connection (with male or female thread)

- Manually screw the thread onto the cylinder valve connector.
- Turn the regulator into the correct position for use and tighten the nut by hand - do not use tools.
- Position the equipment in such way so that the regulator user outlets point away from personnel.

-  **Fitting the regulator to the cylinder valve with too high torque may result in damage.**

-  **When fitting the regulator to the cylinder valve, do not apply torque/ load to any other parts of the product.**

Leakage check before use

- Slowly open the cylinder valve by turning the hand wheel in anticlockwise direction approx. 1 to 1½ turns.

-  **Sudden opening of the cylinder valve could cause danger of fire or explosion arising from oxygen pressure shock. Insufficient opening of the cylinder valve could reduce actual flow delivered.**

- Check visually and acoustically for possible leakages on:
 - the regulator inlet connection to the cylinder valve,
 - the pressure relief valve vent hole(s),
 - the pressure gauge connection to the main body,
- Turn o the cylinder valve by turning the hand wheel in a clockwise direction to the "stop" position. Do not use excessive force.

-  **If any leakage is detected, use the procedure in chapter 6.3 and return the product to GCE for service.**

Functional checks before use

- Ensure that the cylinder valve is open, i.e. in the "ON" position.
- Check that the gauge indicates pressure.
- Turn off the cylinder valve by turning the hand wheel in a clockwise direction to the "stop" position. Do not use excessive force.
- For regulators fitted with a pressure outlet, ensure it is functional by connecting and disconnecting an output extension.

6.2. USER OUTLET(S) CONNECTION & USE

List of recognised accessories

Hoses, flow meters, ventilators, low pressure regulators.

To be connected to the pressure outlet:

Ventilators

-  **Before connecting any accessory or medical device to the regulator, always check that it is fully compatible with the product connection features & the product version.**

EN

Pressure outlet connection

Pressure outlet

- Ensure that the extension is compatible with the pressure outlet.
- Connect the extension.
- Check correct connection of the extension.

-  **Regulators with threaded connector as pressure outlet can only be used as an integral part of a medical equipment unit. Do not use them for other purposes!**

-  **Sudden opening of the cylinder valve could cause danger of fire or explosion arising from oxygen pressure shock. Insufficient opening of the cylinder valve could reduce actual flow delivered.**

-  **The oxygen flow rate must be prescribed and administered by a clinically trained user.**

Use of the product pressure outlet

- Ensure that the flow control knob is in the ZERO position (only for products with flow measuring devices).
- Ensure that the accessory IS NOT connected to the pressure outlet.
- Slowly open the cylinder valve by turning the hand wheel in anticlockwise direction approx. 1 to 1½ turns.

-  **Sudden opening of the cylinder valve could cause danger of fire or explosion arising from oxygen pressure shock. Insufficient opening of the cylinder valve could reduce actual flow delivered.**

6.3. AFTER USE

- Turn off the cylinder valve by turning the hand wheel in a clockwise direction to the „stop“ position. Do not use excessive force.
- Ensure that the flow control knob is set to the “ZERO” value until the flow control device reaches the correct position (only for products with flow measuring devices).
- Check that the pressure gauge does not indicate any residual pressure.
- Disconnect all devices connected to the user outlets.
- Refit pressure outlet and flow outlet protection caps. Before refitting the caps, ensure that they are clean.

7. CLEANING

Remove general contamination with a soft cloth damped in oil free, oxygen compatible, soapy water, & rinse with clean water.

Disinfection can be carried out with an alcohol-based solution (spray or wipes). If other cleaning solutions are used, check that they are not abrasive and that they are compatible with the product materials (including labels) and gas.

-  **Do not use cleaning solutions containing ammonia!**
-  **Do not expose the device to influences of water or any other liquids.**
-  **Do not expose to high temperatures (such as in an autoclave).**

8. MAINTENANCE

8.1. SERVICE AND PRODUCT LIFE TIME

8.1.1. Service

GCE recommend that a Periodic inspection of the product is undertaken every year, including check of proper functionality of the regulator. It should be done by a skilled technician.

GCE recommend that Overall maintenance is performed after 5 years of operation. Such maintenance consists of preventive operational maintenance including replacement of critical components and re-testing of the product. Overall maintenance can be carried out by GCE authorised person only.

We have to warn that the Periodic Inspections and Overall Maintenance recommended by GCE do not necessarily cover every safety procedure or practice required by local regulations or statutory requirements and that abnormal or unusual circumstances can cause further requirements or additional procedures.

8.1.2. Life time and waste management

Maximum life time of the product is 10 years.

At the end of the product's life time, the product must be taken out of service. The device owner must prevent any future use of the product (marking, etc.).

The provider of the device shall prevent the reuse of the product and handle the product in compliance with "Directive of European Parliament and Council 2008/98/EC on waste".

In accordance to Article 33 of REACH GCE, s.r.o. as responsible manufacturer shall inform all customers if materials containing 0.1% or more of substances included in the list of Substance of Very High Concern (SVHC).

The most commonly used brass alloys used for bodies and other brass components contain 2-3% of lead (Pb), EC no. 231-100-4, CAS no. 7439-92-1. The lead will not be released to the gas or surrounding environment during normal use. After end of life the product shall be scrapped by an authorized metal recycler to ensure efficient material handling with minimal impact to environment and health.

To date we have no information that indicates that other materials containing SVHC of concentrations exceeding 0.1% are included in any GCE product.

8.2. REPAIRS

Repair activities cover the replacement of the following damaged or missing components:

- Inlet stem
- Diaphragms of the 1st a 2nd level
- Inlet pressure gauge
- Throttling devices of the 1st a 2nd level
- Outlet pressure gauge
- Pressure relief valve
- Outlet connector

Repairs can be carried out by a GCE authorised person only.

Any product sent back to a GCE authorised person for maintenance must be properly packaged.

The reason for the maintenance has to be clearly specified (repair, overall maintenance). A short description of the fault and the claim number reference should accompany any product sent for repair.

Some repair activities concerning replacement of damaged or missing components can be carried out by the owner of the product. Only the following parts can be replaced:

- Caps
- Hose nipple
- Stickers
- Inlet stem seals

 **All labels on the product must be kept in good, legible condition by the owner and the user during the entire product life time.**

 **All seals and o-rings must be kept in dry, dark, and clean environment by the owner and the user during the entire product life time.**

 **Use only original GCE components!**

9. GLOSSARY

EN

	Consult instruction for use		Suitable for Hospital care use
	Humidity limit		Suitable for Home care use
	Caution		Suitable for Emergency care use
	Keep away from heat and flammable material		Temperature limit
	Keep away from oil and grease		Serial number
	Keep dry		Catalogue number
	Fragile, handle with care		Batch code
	Date of manufacture		Manufacturer
	Use by date		Outlet parameter
	Inlet parameter (P1)		

10. WARRANTY

The Standard Warranty period is two years from date of receipt by the GCE Customer (or if this is not known 2 years from time of the product manufacture shown on the product).

The standard warranty is only valid for products handled according to Instruction for use (IFU) and general industry good practice and standards.

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© GCE, s.r.o.